

APPLICATION

The QD Head coil (MJQH-142A/S1) is designed to produce high-SNR images with enhanced patient comfort for Vantage Titan 3T MRI systems.

COMPOSITION


- QD Head coil..... 1 set
 - Head base
 - Anterior piece
- Mirror..... 1
- Mats and pads 1 set
- Operation manual 1

PERFORMANCE SPECIFICATIONS

This coil is a transmit-receive coil with a 16-element bird-cage design. The open design of the coil and the detachable mirror improve patient comfort. The internal diameter of the coil is 27 cm.

APPLICABLE CODES

IEC60601-1 (1988) + Amendment 1 (1991)
+ Amendment 2 (1995)
IEC60601-1-2 (2001) + Amendment 1 (2004)
IEC60601-2-33 (2002) + Amendment 1 (2005)
+ Amendment 2 (2007)




IEC CLASSIFICATION

Class II
Type BF Applied Part

DIMENSIONS AND MASS

Unit	Dimensions W × D × H mm	Mass kg (Approx.)
QD Head coil	463 × 316 × 355	4.8
Others		1.0

OUTLINE DRAWINGS



Unit: mm



TOSHIBA MEDICAL SYSTEMS CORPORATION

<http://www.toshibamedicalsystems.com>

©Toshiba Medical Systems Corporation 2009 all rights reserved.
Design and specifications subject to change without notice.
2009-11 TME/KI

Toshiba Medical Systems Corporation meets internationally recognized standards for Quality Management System ISO 9001, ISO 13485.

Toshiba Medical Systems Corporation Nasu Operations meets the Environmental Management System standard, ISO 14001.

Made for Life and Vantage Titan are trademarks of Toshiba Medical Systems Corporation.