



MAGNETOM ESSENZA

Tim [25x8]
V-engine

www.siemens.com/medical

SIEMENS
medical



Tim – Total imaging matrix

[25x8]. 25 seamlessly integrated coil elements.
8 RF channels.

- Up to 25 simultaneously connected coil elements which can be seamlessly integrated into one examination
- 8 independent receiver channels (Analog/Digital Converters, ADCs)

Almost all receiving coils with up to 25 coil elements in total can be connected simultaneously. They can be seamlessly integrated into the examination without repositioning the patient or even changing a single coil.

iPAT – Advanced iPAT capabilities

- Full iPAT throughout the whole body – without patient repositioning or changing the coil setup
- Multi-directional, i.e., three dimensional, high-speed, high-resolution iPAT
- Flexible iPAT through use of multiple coils and Matrix Coil Mode
- iPAT with acceleration factors up to 4 (one direction) or 8 (with iPAT², optional)¹⁾
- Tim Assistant ensures ease-of-use and optimized iPAT settings

1) With iPAT Extension Option

Gradients

General

Gradient duty cycle	100%
Gradient performance	
Max. amplitude	30 mT/m
Min. rise time	300 µs
Max. slew rate	100 T/m/s
Vector gradient performance (vector summation of all 3 gradient axes)	
Max. eff. amplitude	52 mT/m
Max. eff. slew rate	173 T/m/s

Gradient Amplifier

Water-cooled, highly compact, modular design
Ultra-fast solid-state technology with very low switching losses

Max. output voltage²⁾	1200 V
Max. output current²⁾	150 A

2) Values for each of the 3 gradient axes

Resolution Parameters			
	Min. FoV	5 mm	
	Max. FoV ¹	45 cm	
	Min. slice thickness 2D	0.1 mm	
	Max. slice thickness 2D	200 mm	
	Min. partition thickness 3D	0.05 mm	
	Max. partition thickness 3D	20 mm	
	Min. slab thickness 3D	5.12 mm	
	Max. slab thickness 3D	450 mm	
	Max. matrix	1024	
	Highest in-plane resolution	16 µm	
1) Depending on the application, the maximum FoV may be smaller (e.g., up to 35 cm in the z-direction)			
Spin Echo	Matrix	128	256
	Min. TR	9.7 ms	11 ms
	Min. TE	3.6 ms	4.4 ms
2D GRE (TurboFLASH)	Matrix	128	256
	Min. TR	1.46 ms	2 ms
	Min. TE	0.62 ms	0.89 ms
	Min. measurement time	34 ms	43 ms
3D GRE (ceMRA)	Matrix	128	256
	Min. TR	1.46 ms	1.99 ms
	Min. TE	0.62 ms	0.9 ms
TrueFISP	Matrix	128	256
	Min. TR	2.24 ms	2.96 ms
	Min. TE	0.94 ms	1.25 ms
	Min. measurement time	72 ms	100 ms
TSE (HASTE)	Matrix	128	256
	Min. echo spacing	2.66 ms	3.0 ms
	Min. TR	13 ms	16 ms
	Min. TE	3.7 ms	4.4 ms
	Min. measurement time	71 ms	77 ms
	Max. Turbo factor	512	

EPI (single-shot and multi-shot)	Matrix	64	128	256
	Min. echo spacing	0.4 ms	0.6 ms	0.97 ms
	Min. TR	10 ms	10 ms	10 ms
	Min. TE	2.4 ms	2.6 ms	3.2 ms
	Min. measurement time	10 ms	13 ms	16 ms
	Max. EPI factor	256		
Diffusion Imaging	Matrix	64	128	256
	Max. b-value [s/mm ²]	10 000	10 000	10 000
	Min. TE with b=1 000 [s/mm ²]	66 ms	70 ms	76 ms

All matrices without interpolation

Combinations of the stated parameters are not always possible; some parameters may require optional application packages

© 02.2008, Siemens AG
Order No.
A91MR-410-12-7600
Printed in Germany
MC.MR-00410 PA 02082.

Headquarter
Siemens AG, Medical Solutions
Henkestr. 127, D-91052 Erlangen
Germany
Telephone: +49 9131 84-0
www.siemens.com/medical

Contact in the USA
Siemens Medical Solutions USA, Inc.
51 Valley Stream Parkway
Malvern, PA 19355
Telephone: +1 888-826-9702
Telephone: +1 610-448-4500
Telefax: +1 610-448-2254

www.siemens.com/medical