

# OASIS™

## High-Field, Bore-Less MR



Comfortably meeting  
your MRI challenges

**HITACHI**  
Inspire the Next

# High-Field Imaging Excellence

*"The Hitachi Oasis has essentially made our Open MRI services limitless. Patients are able to tolerate its wide open design very well. The image quality is exceptional. Our radiologists have expressed that, in some cases, the image quality exceeds even that of our 1.5T conventional MRI system."*

Sandra Holman  
Administrative Director  
Frye Regional Medical Center

*"In the first seven months working with the scanner, we have seen a 19% incremental increase in patient volume. The OASIS now handles the bread-and-butter outpatient exams and it can also handle patients who are not candidates for either of the closed-bore 1.5T or 3T magnets at SLU Hospital."*

Jeffrey Dossett  
Director of Imaging Services  
St. Louis University Hospital





## Neuro – Brain

Vital pulse sequences, acquisition features and post processing tools for high-quality imaging of the brain.



## Neuro – Spine

Oasis' standard CTL coil sensitivity and uniformity complements spine imaging sequences and tools.




## MSK – Upper

High SNR potential from Oasis' vertical field and iso-center positioning promotes high spatial resolution critical for orthopedic imaging.




## MSK – Lower

The ability to move the patient table laterally and perform imaging on iso-center gives Oasis inherent advantages in comfort and image quality.



## Body

High SNR from the 1.2T magnet and Zenith™ RF coil technology is complemented by 2D and 3D protocols for abdomen, pelvis, MRCP and liver.



## Breast

When coupled with the 8-channel Breast coil, Oasis' suite of Breast imaging features delivers excellent image quality and broad capability.




## Vascular

Conventional 2D/3D TOF and advanced acquisition techniques such as Time Resolved MRA (TRAQ™) and 3D vessel post-processing are included.



## Cardiac

Basic cardiac imaging is supported by standard dark blood and bright blood sequences and the 6-channel Body coil.



## Pediatric

With its wide open design and RADAR™ motion compensation, Oasis provides the ideal platform for comfortable pediatric imaging.



## Bariatric

Expansive room for the patient, a 660 lb.- capable table, and the largest flex coil in the industry accommodate bariatric cases other systems can't handle.




- RADAR minimizes motion and flow effects.  
Complete brain studies can be done with RADAR
- RAPID™ parallel imaging technique, enables short scan times
- High performance gradients drive fast, high resolution studies
- HOAST™ promotes fat saturation uniformity for small and large FOV
- Diffusion Weighted SS-EPI with seamless ADC map and isotropic DW image creation
- Zenith RF coils drive high resolution imaging with high SNR



# Neuro - Brain



RADAR assures a diagnostic result even when there is extreme patient motion



RAPID and Driven Equilibrium promote fast scans with high SNR



Resolution: <0.6mm  
Slice: 1.0mm  
Scan time: 3:40




Resolution: <1mm  
Slice: 3.0mm  
Scan time: 3:45

Thin slices and high resolution with short scan times



HOAST delivers excellent fat saturation



DWI ADC

RAPID reduces susceptibility artifacts on SS-EPI diffusion



Large field of view, high resolution and high signal uniformity with Oasis NeuroVascular Coil





- RADAR minimizes motion and flow effects. Complete spine studies can be done with RADAR
- RAPID parallel imaging technique enables short scan times
- High performance gradients drive fast, high resolution studies
- primeFSE reduces susceptibility artifacts from prostheses
- HOAST promotes fat saturation uniformity for small and large FOV
- BASG delivers excellent nerve root delineation
- Zenith RF coils drive high resolution imaging





Spin Echo - T1W



T1 FLAIR



DE-FSE T2

Complete study can be done using RADAR to minimize artifacts from voluntary or involuntary motion



Driven Equilibrium keeps scan time low for T2 weighting



Fast gradients deliver T1 weighting in a short scan time



HOAST promotes excellent large field of view fat sat saturation



primeFSE minimizes susceptibility artifacts from prostheses




Myelogram at 400ms TE and 80ETL



Nerve root detail with BASG sequence

- Isocenter positioning enabled with lateral movement and wide patient table
- Dedicated RAPID orthopedic imaging coils for high SNR and signal uniformity
- HOAST and regional shimming deliver outstanding extremity RF fat saturation
- RADAR motion compensation
- Driven Equilibrium FSE enables heavy T2 weighting with fast scan time







Oasis' extra-wide patient table can move laterally to position anatomy of interest on iso-center. The result is high image quality with the patient maintaining a comfortable position.



Scan time: 3:33  
Resolution: <0.53mm



Scan time: 5:02  
Resolution: <0.40mm




Scan time: 3:35  
Resolution: <0.40mm

Oasis' powerful gradients promote high resolution imaging with fast scan times



Excellent cartilage depiction with RSSG or BASG and water excitation



Driven Equilibrium for PD weighting with fat saturation

- Isocenter positioning enabled with lateral movement and wide patient table
- Dedicated RAPID orthopedic imaging coils for high SNR and signal uniformity
- HOAST and regional shimming deliver outstanding extremity RF fat saturation
- Water Excitation with 3D Gradient Echo sequences for cartilage imaging
- primeFSE reduces susceptibility artifacts from prostheses
- RADAR motion compensation
- Driven Equilibrium FSE enables heavy T2 weighting with fast scan time





# MSK - Lower



Oasis' extra-wide patient table can move laterally to position anatomy of interest on iso-center. The result is high image quality with the patient maintaining a comfortable position.



Oasis' powerful gradients promote high resolution imaging with fast scan times

primeFSE provides for imaging in the presence of prostheses



Excellent cartilage depiction with RSSG and water excitation

Driven Equilibrium fast IR (FIR) with fast scan time



- TIGRE™ fat suppressed 3D dynamic imaging for thin slice and complete organ coverage
- RAPID reduces breath hold and study times. RAPID 3D accelerates volume acquisitions for dynamic imaging
- RADAR for comprehensive free breathing abdominal series
- HOAST optimizes field uniformity for excellent image quality and fat suppression
- In-phase/out-of-phase results in one breath hold
- Large to small FOV capability promoted by sensitive RF Coils





Pre



Post




Post + 5 min.

TIGRE Gradient Echo Sequence for T1 weighting with fat saturation for dynamic liver maging. Fast, 18-second breath hold



Excellent uniformity over large fields of view



MRCP with heavy T2 weighting and MIP processing



Free breathing kidneys with RADAR



Dual echo gradient echo provides in-phase and out-of-phase images in a single, 23-second breath hold



Excellent soft tissue depiction with water excitation



- TIGRE for dynamic studies meets ACR guidelines
- RAPID 3D further accelerates volume acquisitions
- Sensitive Breast RF Coil delivers bilateral, high SNR coverage and patient comfort
- Workflow-focused image processing
- Automatic position correction and subtraction
- DICOM exportable time/intensity curves
- Compatibility with CAD







# Women's Health



Driven Equilibrium T2 Weighting




Maximum Intensity Projection



Pre

8-channel bi-lateral breast coil provides for high SNR and RAPID parallel imaging



Dynamic Tissue Intensity analysis. Images can be sent by DICOM to a separate CAD workstation



Post



Post + 1 minute



RADAR reduces motion artifacts



Uniform fat saturation with HOAST


High resolution pelvic imaging

Dynamic acquisition in full compliance with ACR Guidelines:

- Sub-millimeter in-plane resolution
- <3mm slice thickness
- Bilateral imaging
- Repeat series in 3 min or less (typically <90 sec)

- 2D and 3D TOF for intra-cranial and extra-cranial imaging
- FLUTE™ ensures arterial phase capture
- TPEAKS™ centric k-space filling
- TRAQ MRA depicts blood flow dynamics
- PVA with multi-coil connections, auto table movement, and large FOV
- VASC™ non-contrast MRA for renal, carotid and peripheral applications
- Dedicated RF Coils provide outstanding coverage and uniformity
- Volume Rendering and auto MIP preview






FLUTE, fluoro triggering, and TPEAKS ensure excellent depiction of the arterial phase



VASC for high quality MRA without a contrast agent

Volume rendering helps visualize vascular abnormalities

Auto Table Step enables multi-station runoff studies. Image stitching creates a single image for archive or export



Sensitive Zenith RF coils provide for high SNR and uniformity and large field of view



Sensitive Zenith coils provide detail in lower vasculature



- Flexible, workflow-focused tools to assess cardiac morphology and function
- Double and triple inversion recovery breath hold acquisitions
- Navigator echo free-breathing scan
- RAPID parallel imaging increases temporal resolution and image quality
- Multi-slice, multi-phase cine
- Viability assessment benefits from powerful gradients and sensitive RF Coils





Double Inversion Recovery





Triple Inversion Recovery






Free breathing navigator echo



Multi-slice/multi-phase cine for functional assessment



Myocardium assessment with delayed enhancement




Viability assessment with dynamic tissue intensity



- Imaging free of ionizing radiation – a safer alternative for pediatric patients
- Constant parent-child contact dramatically improves patient compliance
- RADAR motion free imaging eliminates repeat scans and improves image quality
- Fast scanning techniques keep study time minimal
- Halo coil delivers quality imaging and an all-around view









Oasis' RADAR feature delivers diagnostic results even with non-compliant patients



White matter suppression with primeFIR



High quality pediatric head imaging with the solenoid "halo" coil



Echo Planar sequence for DW weighting with ADC



FOV: 14cm  
Resolution: <0.73mm

Small FOV – High resolution MSK imaging



Out of phase Gradient Echo



Aortic arch

High SNR Body Imaging

- Oasis' coil offerings accommodate larger patients while maintaining high SNR
- NATURAL™ provides signal intensity compensation necessary for larger patients
- Optimal image quality with iso-center positioning for even the largest of patients
- 82 cm wide patient table with 660 lb. capacity
- Extra large flexible body coil for patients at the extreme of the demographic spectrum
- Integrated transmit/receive coil for good imaging results if no other coil is appropriate








NATURAL provides signal compensation necessary in spine studies for larger patients




Scan Time: 3:43

Driven Equilibrium combines with prime-FIR to keep scan time low



FOV: 40cm

High quality body imaging with the Extra Large flex body coil



MSK imaging of large patients is less challenging with Oasis due to the wide open design, large table and iso-center positioning




# Image Annotation Guide

Locations and meaning of annotation text used with Oasis images

Patient Name  
Sex/Age/Weight  
Slice#/Total Slices

Site Name  
Date/Time

Contrast



Phase Encoding Direction  
Scan Time  
TR/TE  
Flip Angle  
TI Time  
FOV/Slice Thickness/Interval  
Freq x Phase/RAPID/NSA/Recon Matrix  
Seq. Type/Shot #/Echo Factor/Fat Sat  
b-factor

Sequence Name  
Slice Position

Window  
Level  
System Name

# OASIS™

**Providing the images  
you need and  
the comfort your  
patients deserve**



**Critical Care**



**Pediatric**



**Anxious**



**Open Access**



**Bariatric**

# **HITACHI**

## **Inspire the Next**

### **Hitachi Medical Systems America, Inc.**

1959 Summit Commerce Park  
Twinsburg, Ohio 44087 USA  
Tel: 330.425.1313 800.800.3106  
Fax: 330.425.1410  
[www.hitachimed.com](http://www.hitachimed.com)

### **Hitachi Medical Corporation**

4-14-1 Akihabara UDX  
Soto-Kanda, Chiyoda-ku  
Tokyo, 101-0021 Japan  
[www.hitachi-medical.co.jp](http://www.hitachi-medical.co.jp)

©Hitachi Medical Systems America, Inc. 2010.  
All rights reserved.  
0510/1000/DM# 71403 v1