



MAGNETOM Sola

Tim [204x48] XQ Gradients

Magnet System

General

Superconducting Magnet	Short bore, patient-friendly design, high homogeneity 1.5 Tesla with 70 cm Open Bore design
	Easy siting due to AS (Active Shielding) and E.I.S. (External Interference Shielding) magnet technology

Magnet Parameters

Operating field strength	1.5 Tesla
Magnet type	Superconductor
Field stability over time	< 0.1 ppm/h
Weight (with cryogenics)	2700 kg
Magnet length	1.45 m
Open Bore design¹⁾	70 cm
System length cover to cover	1.57 m

Homogeneity (based on highly accurate 24 plane plot)

10 cm DSV	Guaranteed	0.006 ppm
	Typical	0.004 ppm
20 cm DSV	Guaranteed	0.05 ppm
	Typical	0.03 ppm
30 cm DSV	Guaranteed	0.20 ppm
	Typical	0.11 ppm
40 cm DSV	Guaranteed	0.75 ppm
	Typical	0.65 ppm
50 cm DSV	Typical	5.5 ppm
50 × 50 × 45 cm² DEV	Guaranteed	3.0 ppm
	Typical	2.8 ppm

In compliance with the German "Qualifikationsvereinbarung".

Standard deviation Vrms (volume root-mean square) measured with highly accurate 24 plane plot method (20 points per plane).

Standard active shim with 3 linear channels.

DSV = Diameter spherical volume; DEV = Diameter elliptical volume (x, y, and z directions).

¹⁾Incl. shim coils, gradient coil, RF body coil

Magnet System

Shimming		
Both: passive and active shimming. Passive shimming during installation		
Standard active shim with 3 linear channels (1 st order) and 5 non linear channels (2 nd order)		
3D Shim	Patient-specific automated shim	
	Time to shim	Approx. 15 s
Shielding		
Active Shielding (AS)	5 th generation active shielding (AS) technology with counter coils	
Fringe field (axial × radial)	0.5 mT ¹⁾	4.00×2.50 m
	0.1 mT	5.8×3.4 m
External Interference Shield (E.I.S.)	Patented shielding system integrated into the magnet	
	Continuous compensation and automatic suppression of external magnetic field interferences during measurement (caused by moving ferromagnetic objects or nearby power lines)	
Magnet Cooling System		
Zero Helium boil-off technology		
Refill interval (typical) ²⁾	Not applicable	
Boil-off rate (typical) ²⁾	0.0 liter/year	

¹⁾ Pacemaker safety limit

²⁾ For typical clinical use, depending on pulse sequences and operating time with running helium compressor. The system needs to be serviced at regular interval. Undisturbed magnet cooling for 24 hours and 7 days a week.

XQ Gradients

General Features

- Actively shielded (AS) whole-body gradient coil system
- Extremely low eddy currents
- Water-cooled coil and amplifier for maximum performance
- All axes force compensated

Gradient Performance for Each Axis

Max. amplitude 45 mT/m

Min. rise time 225 μ s

Max. slew rate 200 T/m/s

Vector Gradient Performance (vector addition of all 3 gradient axes)

Max. eff. amplitude 78 mT/m

Max. eff. slew rate 346 T/m/s

Gradient duty cycle 100%

Resolution Parameters

Min. FoV 5 mm

Max. FoV¹⁾ 500 mm

Slice thickness 2D min. 0.1 mm, max. 200 mm

Partition thickness 3D min. 0.05 mm, max. 20 mm

Slab thickness 3D min. 5 mm, max. 500 mm

Max. matrix 1024

Highest in-plane resolution 12 μ m

¹⁾ Depending on the application, the maximum FoV in the z-direction can be up to 50 cm

XQ Gradients

Power Consumption ¹⁾	
System Off	4.3 kW
System ready to measure ¹⁾	8.7 kW
Scan ¹⁾	22.7 kW
Gradient Amplifier	
Water-cooled, highly compact, modular design	
Ultra-fast solid-state technology with very low switching losses	
Max. output voltage ²⁾	2250 V
Max. output current ²⁾	900 A
Max. power ²⁾	2.025 MVA
Power Requirements	
Line Voltage	380V, 400V, 420V, 440V, 460V, 480V
Stability tolerance	± 10 %
Line frequency	50/60 Hz, ± 1Hz
Connection Value	88 kVA

Cooling

Two different customer specific cooling alternatives (Separator or Eco Chiller) available

Eco Chiller option with automatic adaptation to the required cooling demands (e.g. different night/day mode) to decrease energy cost

GREEN Cooling Package³⁾: The Free Cooling Unit reduces energy consumption by up to 39%. It automatically starts if the surrounding temperature is 18°C (64°F) or less and reduces the chiller energy consumption. If the temperature is less than -10 °C (14°F), the chiller is switched off⁴⁾.

Cooling with separator

Water consumption	100 l/min +/- 10 l/min ⁵⁾
Heat dissipation to water	60 kW

¹⁾ All values are typical values, applicable for 400V/50Hz.

The power consumption measurement is based on the COCIR methodology – MRI – Measurement of energy consumption. Many variables impact power consumption, thus there can be no guarantee that each customer will achieve the same values.

²⁾ Values for each of the 3 gradient axes

³⁾ Optional: based on climatic dates of Munich; data on file; results may vary

⁴⁾ In case of clinical routine measurement conditions

⁵⁾ Water temperature 6°-14°C; allowed delta T: +/- 2 K with max. 1K/30s

XQ Gradients

Sequences:

		Matrix		
		64	128	256
Spin Echo	min. TR [ms]	5	5	5.5
	min. TE [ms]	1.5	1.5	1.8
Inversion Recovery	min. TR [ms]	26	26	26
	min. TE [ms]	1.5	1.5	1.8
	min. TI [ms]	21	21	21
2D GRE	min. TR [ms]	0.59	0.7	0.97
	min. TE [ms]	0.22	0.22	0.22
3D GRE	min. TR [ms]	0.59	0.7	0.97
	min. TE [ms]	0.22	0.22	0.22
TrueFISP	min. TR [ms]	1.66	1.81	2.28
	min. TE [ms]	0.76	0.77	0.98
TSE (HASTE)	min. Echo Spacing [ms]	1.54	1.54	1.84
	min. TR [ms]	5	5	5.5
	min. TE [ms]	1.5	1.5	1.8
	max. Turbo Factor = 512			
Turbo GSE	min. Echo Spacing [ms]	0.78	0.82	0.86
	min. TR [ms]	5.7	6.3	6.4
	min. TE [ms]	3	3.5	3.5
	max. Turbo Factor = 65			
	max. EPI Factor = 21			
EPI (single-shot and multi-shot)	min. Echo Spacing [ms]	0.28	0.49	0.66
	min. TR [ms]	10	10	10
	min. TE [ms]	2.1	2.3	2.7
	min. Measurement time	11	17	26
	max. EPI Factor = 256			
Diffusion Imaging	Max. b-value [s/mm ²]	10 000	10 000	10 000
	Min. TE [ms] with b = 1000 s/mm ²	40	42	46

All matrices without interpolation. Combinations of the stated parameters are not always possible; some parameters may require optional application packages.

Coils

1.5T BioMatrix Coils

A new ultra-high density BioMatrix coil (BioMatrix Spine 32) utilizes seamlessly integrated sensors to acquire and display the patient's respiration data without need for user interaction.

The integrated CoilShim technology in the BioMatrix Head/Neck 20 ensures that the challenging head/neck region is automatically and optimally shimmed for reproducible quality in every patient.

1.5T Tim 4G Coils

The Tim 4G coils are designed for highest image quality in combination with easy handling. High element coils increase SNR and reduce examination times. DirectConnect® and SlideConnect® technology reduce patient set up time. Light weight, ergonomically designed coils enable highest patient comfort.

- No coil changing with multi-exam studies saves patient setup time
- All coils are time-saving "no-tune" coils
- Low-noise preamplifiers
- AutoCoilSelect for dynamic, automatic, or interactive selection of the coil elements within the Field of View

Standard Coil Package - Tim [204 x48]

BioMatrix Head/Neck 20 tiltable with CoilShim (DirectConnect)	Application area	Head and Neck	
	Dimensions (L × W × H)	425 mm × 370 mm × 385 mm	
	Weight	5.7 kg	
BioMatrix Spine 32 with Respiratory Sensors (DirectConnect)	Application area	Spine	
	Dimensions (L × W × H)	1200 mm × 489 mm × 63 mm	
	Weight	10.5 kg	
Body 18 (SlideConnect)	Application area	<ul style="list-style-type: none">• Thorax• Heart• Abdomen	<ul style="list-style-type: none">• Pelvis• Hip
	Dimensions (L × W × H)	385 mm × 590 mm × 65 mm	
	Weight	1.6 kg	
Flex Large 4	Application area	Multi purpose	
	Dimensions (L × W)	516 mm × 224 mm	
	Weight	550 g	
Flex Small 4	Application area	Multi purpose	
	Dimensions (L × W)	366 mm × 174 mm	
	Weight	450 g	
Accessories	<ul style="list-style-type: none">• Flex Coil Interface 1.5T		

Combination of all coils possible for large Field of View exams.

Computer System

syngo Acquisition Workplace - Tim [204 x 48]

Standard measurement

and reconstruction system¹⁾

Processor	2x Intel Xeon ≥ E5-2609v4 (8 Core)
Clock rate	2x ≥ 1.7 GHz, or comparable
Main memory (RAM)	64 GB
SSD for raw data	≥ 480 GB
SSD for system software	≥ 240 GB
Parallel Scan & Recon	Simultaneous scan and reconstruction of up to 12 data sets
Reconstruction speed	32 854 recons per second (256 ² FFT, full FoV) 122 137 recons per second (256 ² FFT, 25% recFoV)

Highend measurement and reconstruction system¹⁾

Processor	2x Intel Xeon ≥ E5-2620v4 (8 Core)
Clock rate	2x 3.3 GHz, or comparable
Main memory (RAM)	96 GB
SSD for raw data	≥ 480 GB
SSD for system software	≥ 240 GB
Parallel Scan & Recon	Simultaneous scan and reconstruction of up to 12 data sets
Reconstruction speed	40 404 recons per second (256 ² FFT, full FoV) 149 532 recons per second (256 ² FFT, 25% recFoV)
GPGPU	1× Nvidia Quadro P2000

¹⁾ Optional

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products included in this brochure are available through the Siemens sales organization worldwide. Siemens reserves the right to modify the design, packaging, specifications, and options described herein without prior notice.

Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

Please find fitting accessories:

healthcare.siemens.de/accessories-oem-equipment

Please contact your local Siemens sales representative for the most current information.

International version.

Not for distribution in the U.S.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84 0
siemens.com/healthineers

Legal Manufacturer

Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany